

# Direct access to mass spectrometry data in R - The MsBackendRawFileReader

Tobias Kockmann, Christian Panse

Functional Genomics Center Zurich, Swiss Federal Institute of Technology Zurich | University of Zurich, Winterthurerstr. 190, CH-8057 Zurich, SWITZERLAND.

## Abstract

The R for Mass Spectrometry (RforMS) initiative aims to provide efficient, thoroughly documented, tested and flexible software for the analysis of mass spectrometry data in R (<https://www.rformassspectrometry.org>). Reaching this goal at first requires performant and genuine access to the data from within the R environment. Unfortunately, proprietary vendor formats have so far obstructed this direct access. Our R package **rawR** (source: <https://github.com/cpanse/MsBackendRawFileReader>) eliminates this limitation for Thermo Fisher Scientific raw files by wrapping methods defined in the RawFileReader.NET assembly [1]. The resulting on-disc MsBackendRawFileReader is fully compliant with the **Spectra** class as defined by the RforMS initiative and allows highly efficient and straightforward data mining directly in R.

## 1 Test data

The Bioconductor **tartare** package [2] provides raw files (spectra) recorded on recent Thermo Scientific mass spectrometers.

```
R> library(tartare)
R> eh <- ExperimentHub()
R> query(eh, c('tartare'))

ExperimentHub with 4 records
# snapshotDate(): 2019-10-22
# dataProvider: Functional Genomics Center Zurich (FGCZ)
# $species: NA
# $dataType: Spectra
# additional $cols(): taxonomyid, genome, description, coordinate_1_based, maintainer,
# $rdatadate, preparerclass, tags, rdathpath, sourceuri, sourcetype
# retrieve records with, e.g., 'object[["EH3219"]]'

title
EH3219 | Q Exactive HF-X mzXML
EH3220 | Q Exactive HF-X raw
EH3221 | Fusion Lumos mzXML
EH3222 | Fusion Lumos raw
```

```
R> files <- getFilename(eh)
```

## 2 Implementation

The **rawR** codebase operates on multiple language and hierarchy levels. High level R functions call C# definitions to directly access mass spectrometry data stored in vendor-formatted binary files (\*.raw). These levels are connected by an in-process bridge (rDotNet) for the .NET framework (see Figure 1). This section gives an overview of the **rawR** C# definitions and their execution from the R command line.




Figure 1: The graph depicts the embedding into the RforMS packages. The node [.NET] is explained in section 2.1. The proof-of-concept of the MsBackendRawFileReader is drafted in section 2.2. The ".NET" libraries are loaded by a runner CLRServer.exe that provides a [TCP-based API], giving full visibility into your library(es) [3, README]. An application on using the Spectra class with rawDiag [4] is demonstrated at the end.

### 2.1 Calling the vendor library from C#

The following listing demonstrates the interaction with the vendor specific libraries (ThermoFisher.CommonCore\*.dll) in C#. The vendor libraries provide access to the raw file structure.

```
12 using System;
13 using System.Collections.Generic;
14 using System.Diagnostics;
15 using System.IO;
16 using System.Runtime.ExceptionServices;
17 using System.Collections;
18 using System.Linq;
19 using ThermoFisher.CommonCore.Data;
20 using ThermoFisher.CommonCore.Data.Business;
21 using ThermoFisher.CommonCore.Data.FilterEnums;
22 using ThermoFisher.CommonCore.Data.Interfaces;
23 using ThermoFisher.CommonCore.MassPrecisionEstimator;
24 using ThermoFisher.CommonCore.RawFileReader;

27 namespace MsBackendRawFileReader
28 {
29     public class Rawfile
30     {
31         private IRawDataPlus rawFile;
32         public Rawfile(string rawfile)
33         {
34             // ...
35         }
36     }

127     public bool IsValidFilter(string filter)
128     {
129         // ...
130     }

137     public string[] GetAutoFilters()
138     {
139         // ...
140     }

398     public double GetRTinSeconds(int scanNumber) {
399         // ...
400     }

551     public double[] GetSpectrumNoises(int scanNumber, string scanFilter)
552     {
553         // ...
554     }
555 }
```

At the moment there are 40 "getter" methods implemented. Below we list all method signatures used on this poster:

```
563     public double[] GetSpectrumCharges(int scanNumber, string scanFilter)
564     {
565         // ...
566     }

589     public double[] GetSpectrumResolutions(int scanNumber, string scanFilter)
590     {
591         // ...
592     }

603     public double[] GetSpectrumIntensities(int scanNumber, string scanFilter)
604     {
605         // ...
606     }

620     public double[] GetSpectrumMasses(int scanNumber)
621     {
622         // ...
623     }
624 }
```

### 2.2 Calling the C# methods from R

The snippet below calls the constructor for a given rawfile.

```
R> library(MsBackendRawFileReader)
R> ## https://CRAN.R-project.org/package=rDotNet
R> ## create an object and call a method
R> rawfile <- files[4]
R> (x <- rDotNet::cnew("Rawfile", rawfile))

<dotnet obj: 4, class: MsBackendRawFileReader.Rawfile, value: "MsBackendRawFileReader.Rawfile">
```

An implementation in C# of extracting the mass over charge (m/z) values of a spectrum is listed below.

```
630     public double[] GetSpectrumMasses(int scanNumber)
631     {
632         var scanStatistics = rawFile.GetScanStatsForScanNumber(scanNumber);
633         var centroidStream = rawFile.GetCentroidStream(scanNumber, false);
634         if (scanStatistics.IsCentroidScan && centroidStream.Length > 0)
635             return centroidStream.Masses.ToArray();
636         else
637         {
638             var segmentedScan = rawFile.GetSegmentedScanFromScanNumber(scanNumber, scanStatistics);
639             return segmentedScan.Positions.ToArray();
640         }
641     }
642 }
```

Once the .NET library is registered their methods can be called from R as follows: "The R api uses the \$ syntax to reference members much like other R object approaches." [3, README]. Here is a code snippet in R extracting peak attributes from a single spectrum.

```
R> # Define scan number 4034
R> scan <- 4034
R> # getting a vector of mass values
R> head(mz <- x$GetSpectrumMasses(scan))

[1] 110.0709 111.0744 114.2004 120.0805 129.1020 138.0659

R> # getting a vector of intensities
R> head(intensity <- x$GetSpectrumIntensities(scan))

[1] 328132.188 11543.382 4682.475 20583.205 7574.656 6460.671

R> # getting data not available from mzXML
R> noise <- x$GetSpectrumNoises(scan)
R> resolution <- x$GetSpectrumResolutions(scan)
R> charge <- x$GetSpectrumCharges(scan)
```

```
R> # https://CRAN.R-project.org/package=protViz
R> pp <- protViz::peakplot("HTFSGVASVSSGGEAFHVGK", list(mz = mz, intensity = intensity),
R> FUN = function(b,y){cbind(b,y)}, stol = 0.01, col = 'lightgrey')
R> legend("top", "", title=x$ScanFilter(scan), cex=1.5, bty = "n")
R> lines(mz, noise, col="FF7F0E", lwd=2)
R> points(mz,intensity, pch=16, col = as.factor(charge), cex=0.5)
R> legend("topleft", paste(unique(charge), "*"), pch=16, col=as.factor(unique(charge)))
```




Figure 2: The graphic depicts the peptide spectrum match (PSM) with the "HTFSGVASVSSGGEAFHVGK" peptide using <https://CRAN.R-project.org/package=protViz> of scan 4034. The peptide assignment was determined by precursor mass and a scoring function based on counting fragment ion matches. The graphics on the bottom plot the resolution and noise characteristics as a function of mass over charge.

Scan collections can be subset using native filters on the vendor library level, see also [1, page 28].

```
R> head(x$GetAutoFilters())

[1] "FTMS + c ESI Full ms [350.0000-1500.0000]"
[2] "FTMS + c ESI Full ms [350.0000-1800.0000]"
[3] "FTMS + c ESI Full ms [350.0000-2000.0000]"
[4] "FTMS + c ESI Full ms [375.0000-1500.0000]"
[5] "FTMS + c ESI d Full ms2 376.69138hc430.00 [110.0000-764.0000]"
[6] "FTMS + c ESI d Full ms2 378.23598hc430.00 [110.0000-767.0000]"

R> x$IsValidFilter("most likely not a valid filter string")

[1] FALSE

R> filterString <- "FTMS + c ESI Full ms [375-1500]"
R> x$IsValidFilter(filterString)

[1] TRUE

R> # determine the peptide weight
R> (mZ <- (protViz::parentIonMass("HTFSGVASVSSGGEAFHVGK") + 1.008) / 2)

[1] 1060.506

R> X <- x$GetXIC(mZ, tol = 5, filter = filterString)

intensity
2500000
1000000
0
20 21 22 23 24
time
```




Figure 3: XIC of mass 1060.5061 Da within a mass window of 5 ppm applying the filter FTMS + c ESI Full ms [375-1500]. The red line marks the retention time of 20.2 seconds when the MS2 (scan 4034) shown in Figure 2 was taken.

## 3 Usage

We decided to use the object-oriented representation of spectral data as defined by the RforMS initiative in the **Spectra** package. The core idea is to *untangle* functionality and data origin/storage (see [https://github.com/jorainer/swemsa\\_2019](https://github.com/jorainer/swemsa_2019)). To demonstrate the power of this concept we cross compare data access through MsBackendMzR with MsBackendRawFileReader.

```
R> BioManager::install('cpanse/MsBackendRawFileReader')

R> # initialize bioconductor tartare package data
R> be <- lapply(files, function(f){
+   if (grepl("mzXML", f))
+     backendInitialize(Spectra::MsBackendMzR(), files = f)
+   else
+     backendInitialize(MsBackendRawFileReader(), files = f, extra=TRUE)
+ })
R> S <- lapply(be, Spectra)
R> # Howto merge different backends?
R> sapply(S, class)

"EH3219" "EH3220" "EH3221" "EH3222"
"Spectra" "Spectra" "Spectra" "Spectra"
```

The following lines plot the scan frequency on all recorded MS level of all four files using **rawDiag** [4] plot methods on the **Spectra** objects bound to different backends.

```
R> library(rawDiag)
R> gpl <- lapply(S, function(x){ rawDiag::Plot(x, FUN = "ScanFrequency")})
```




Figure 4: Coplots compare scan frequencies extracted from two raw files accessed by two backend implementations (2x2 design). Left-to-right: Q Exactive mzXML and raw; Fusion Lumos mzXML, raw.

## 4 Discussion and outlook

As shown in Figure 5 the MsBackendMzR exposes only a small number of scan metadata and reads those into memory during the backend initialization process.




Figure 5: The table displays the scan meta data of two Orbitrap instruments using raw and mzXML files. The availability of a variable is indicated by a black rectangle.

In contrast, the MsBackendRawFileReader reads only a few "run-level" metadata, e.g., instrument id or file name during the initialization. The different philosophies in the backend implementations are also visible in the benchmark below.




Figure 6: The boxplots show ten repetitions of initializing and extracting the precursor charge information of a batch of 1:17 raw and mzXML files on a single-core computer. The scatter plot on the right compares the memory size in Bytes of the objects returned by both backendInitialize function calls.

In the near future we hope to fully synchronize **rawR** with low-level infrastructure packages like **Spectra** and **Chromatograms**. Once a stable interoperability is achieved, it will be submitted to Bioconductor in order to disseminate **rawR** into the mass spectrometry community.

## References

- Jim Shofstahl. Reading Files using the C# file reader, December 2015. URL: <http://planetorbitrap.com/rawfilereader#.WjkiUtJmL4>.
- Christian Panse and Tobias Kockmann. *tartare: Raw ground spectra recorded on Thermo Fisher Scientific mass spectrometers*, 2019. R package version 1.0.0. URL: <https://bioconductor.org/packages/tartare/>.
- Jonathan Shore. *rDotNet: Low-Level Interface to the '.NET' Virtual Machine Along the Lines of the R C/Call API*, 2017. R package version 0.9.1. URL: <https://CRAN.R-project.org/package=rDotNet>.
- Christian Trachsel, Christian Panse, Tobias Kockmann, Witold E. Wolski, Jonas Grossmann, and Ralph Schlappach. *rawDiag: An R Package Supporting Rational LC-MS Method Optimization for Bottom-up Proteomics. Journal of Proteome Research*, 17(8):2908-2914, July 2018. URL: <https://doi.org/10.1021/acs.jproteome.8b00173>, doi:10.1021/acs.jproteome.8b00173.